

SOMAerect Vacuum Erection Device Prescription Advice form

Hospital:

DELIVERY:

Dear Doctor: Date:

Regarding patient:

Patients are requested to take their prescription to their local community pharmacy who should order directly on

Tel: 01923 23 77 95

or Fax: 01923 804 206
orders@lmedicare.co.uk

Available through wholesalers
(Special department)
FREE delivery to the Pharmacy
from IMEDicare LTD.

Surname:	Hospital No.:
First names(s):	DOB:
Address:	
SPACE FOR PATIENT IDENTIFICATION LABEL	

Your patient was seen in the outpatient department at this Hospital today. Following physical assessment and full instruction in correct usage it has been established that the device / ring options recommended below will prove most suitable and effective for him as a solution to his erectile dysfunction. Can you please prescribe:

ERECTILE DYSFUNCTION MANAGEMENT SYSTEMS

(Vacuum Therapy Devices or Accessories) To GP - Please mark prescription 'SLS' (Selected List Scheme)

NHS Order Code:

5 year Warranty

SOMACorrect Xtra - (PIP code: 317-5593).

(Manually operated one-piece with - 3 sizes (to fit) cylinder options sizing options (L, M and S) to stop scrotal skin and supra-public fat entrapment for comfortable erection generation). Ultra-Ease Ring Loading System - poor manual dexterity. 15 Erection Maintenance Rings: Reuseable: 4 Sizes SureEase (Practice), 3 Sizes Low Tension: Select. 3 Sizes Medium Tension Comfort: Ultimate. 3 Sizes High Tension: Ultra. Two easy removal single use SureRelease Rings (S and E), Instructional DVD and manual, Travel Case. Also Peyronies disease application - IFU requested separately)

15111

SOMAerect Response II XL - (PIP code: 317-5585)

(Manually operated one piece device with 9 rings, extra wide cylinder diameter and loading cone).

15888

SOMAerect Response II (PIP code: 3131182)

(Manually operated one piece device with 9 rings, Large cylinder size and loading cone)

15019

SOMAerect Touch II - (PIP code: 3131190)

(S172.16) (Battery operated one-piece device with 9 rings - Large cylinder size and loading cone, dexterity problems only)

15013

Assist Maintenance System - (PIP code: 3131232)

(Loading cone, applicator collar, 3 Select rings and lubricant).

15109

Ultra Maintenance Ring Set - (PIP code: 317-5619)

(3 sizes High tension - venous leak and narrow penile diameter).

(SIZES: S M L
4 5 6).

15444

High Tension

Ultimate (SureFit) Maintenance Ring Set (PIP code: 317-5601)

(3 sizes medium-high tension - super comfort).

(SIZES: X Y Z)

15222

Select Maintenance Ring Set (PIP code: 3131661).

(3 sizes - low to medium tension)

(SIZES: 4 5 6)

15230

SureEase Maintenance Ring Set - (PIP code: 3131240).

(4 pack - SureEase Comfort - soft practice rings with easy removal grips).

(SIZES: 2 2 3 4)

15085

Low Tension

Yours sincerely

SIGNED.....NAME (Capitals).....

TITLE.....or circle:

Consultant, Associate Specialist, Registrar, Staff Grade, Specialist Nurse Practitioner, Nurse Consultant
Clinical Assistant, SHO.

Directorate: ...UROLOGY / DIABETES / GUM / SEXUAL HEALTH / CARDIOLOGY / OTHER

Schedule 11 Prescription Criteria for Vacuum Therapy - Diabetes - Multiple Sclerosis - Parkinson's Disease - poliomyelitis - prostate cancer - prostatectomy (includes TURP) - radical pelvic surgery - severe pelvic injury - renal failure treated by dialysis or transplant - single gene neurological disease - spinal cord injury - spina bifida.

